

YMW02MA. Windows software for collecting and storing data from moisture analyzers.

YMW02MA is a PC software solution that expands the functions of series MA 45, MA 50C/H and MA 100C/H moisture analyzers.

This program lets you collect data from MA 45, MA 50/MA 100 analyzers and store it in a file for further processing with a spreadsheet program. It also displays drying progress on-line as a "drying curve."

The YMW02MA software also lets you use your PC to configure new drying routines, which you can store on the hard drive and download as needed to the MA 45, MA 50/MA 100.

Conversely, drying programs can also be loaded from the moisture analyzers into the PC; for example, for editing or for downloading to another moisture analyzer.

System requirements:
PC with Windows 95/98/NT or compatible OS

2 MB available space on the hard disk for program installation

At least 1 MB of additional hard disk space for storage of data collected from moisture analyzers and/or user-defined drying programs

Mouse and data-transfer interface

Functions:
Data collection and storage in separate files
Temporary on-line display of the drying progress as a graphic curve

Creation of new drying routines for MA 45, MA 50C/H & MA 100C/H

Uploading of drying programs from MA 45, MA 50/MA 100 and downloading of programs to MA 45, MA 50/MA 100 moisture analyzers

Equipment supplied:
Program CD

Interface cable

DB9/DB25 interface cable

Operating instructions included in program as on-line help

Choice of English or German program version

Sartorius AG
Weender Landstrasse 94-108
37075 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius.com

Publication No.: W--1091-e02071

